

Another Solution for Cuff Tear Arthropathy

One Stem. More Options.



equinox[®]

CTA Head

ANATOMICAL. REDEFINED.

The Equinoxe® CTA Head enables a surgeon to use the platform stem for the entire continuum of care.

KEY FEATURES

- Extended articular surface to articulate with acromion in cuff deficient patient
- Platform stem facilitates conversion to and from a Reverse TSA
- Six diameters in two heights provide intra-operative flexibility
- Soft-tissue sparing AP geometry
- Offset replicator plates enable further lateralization of the humeral head



OVERVIEW TECHNIQUE



1

Prepare humerus and insert stem.



2

Attach Replicator Plate and Torque-Defining Screw.



3

Resect the Tuberosity using the CTA Head Cut Guide placed over the Replicator Plate.



4

Trial the CTA Head to determine appropriate diameter and height.

Refer to the CTA Head Operative Technique (Lit# 718-06-31 Rev. A) for a full description of steps.

Exactech, Inc. is proud to have offices and distributors around the globe. For more information about Exactech products available in your country, please visit www.exac.com.